

Water Heating & Cooling Pumps and Pads

ThermoCare offers a variety temperature controlled water pumps and pads for safe, consistent, heat therapy on the surgery table. Animals can be placed on the pads, which are attached to the heating/cooling circulator pump, during surgery to keep body temperature at the desired level and reduce risk of shock.

ThermoCare offers the Gaymar T-Pumps and Mul-T-Pads, Aquamatic K-Modules and K-Pads, as well as various accessories and adaptors that can accommodate a range pumps and pads. Pads are either reusable and repairable — or disposable — and come in a variety of sizes. Optional slip-on fleece bags are machine washable for easy removal of blood, feces, or other fluids. ThermoCare can also facilitate pump repairs; please call or email info@thermocare.com for more information.

T/Pump and Mul-T-Pads

Sandwich the patient in continuous warmth by connecting two pads to one pump

T/Pump Advanced Features

LED Display Panel and Selectable Temperature Settings

- Clear and intuitive, LED lights indicate current settings
- Classic T/Pump: Select from Cooling, Low Heat or High Heat
- Professional T/Pump: Select from 50°F, 95°F, 100°F or 107°F
- Solid light indicates temperature set point has been reached

Comfortable Handle

- Large grip area for easy transport of unit

Audible and Visual Alarm

- Alerts caregiver to flow, low water or internal issues

On/Standby Indicator

- Yellow light indicates unit has power
- Green light indicates unit is on

Translucent Reservoir

- Determine water level at a glance
- Reservoir has large capacity for less frequent filling
- Fill lines for both warming and cooling (water and ice)

Non-marking Rubber Feet

- Provide stability and shock absorption
- Feet easily realigned to fit T/Pump stand



Extra-Wide Leakproof Cap

- Simply fill with tap water for warming or tap water and ice for cooling and unit is ready for operation. Fill lines are indicated for both water and ice
- Tethered strap prevents cap loss




Secure Hose Attachments

- Permanently attached 10' hoses assure proper, secure connections
- Hoses can not be accidentally lost or discarded



Leakproof Connections

- Patented Clik-Tite® connectors snap into place quickly and easily
- Allows for two pads to be used simultaneously with one unit to treat multiple sites
- Non-restricted water flow



Accessory T/Pump Stand

- Provides ease of transporting to patient's room
- Storage in wire basket for pads and accessories



T/Flash Technology

- Slight temperature change stimulates patient's heat receptors
- Provides the "always there" feeling so the patient knows it's working

Safe, 3 Layer Safety System

- Dual temperature sensors
- Dual microprocessors
- Back-up mechanical thermostat to ensure safe operation

Safe, Effective Therapy

- When used in combination with our versatile Mul-T-Pads®, the T/Pump provides safe, effective localized therapy while easing patient pain and discomfort



| Model | Water Heating Pumps | Price/Case |
|--------|--|------------|
| KM-107 | Aquamatic K-Module 107 Pump for K-Pad w/10' Dual Tubing & Key (See BXK-2 adaptor to connect Gaymar Mul-T-Pads to K-Module Pump) | \$614.00 |
| GTP-K | Gaymar T-Pump 650C w/Colder-style Connectors for K-Pads (See BXK-2 adaptor to use Gaymar Mul-T-Pads) | \$486.00 |
| GTP | Gaymar T-Pump w/Clik-Tite Connectors for 2 pad use w/Mul-T-Pads (See GTA adaptor to use with K-Pads) | \$486.00 |

| Model | Water Heating & Cooling Pumps (Heats & Cools, Tap Water) | Price/Case |
|--------|---|------------|
| TP700C | Gaymar T-Pump w/Colder-Style Connectors for K-Pads (See BXK-2 adaptor to use Gaymar Mul-T-Pads) | \$536.00 |
| TP700 | Gaymar T-Pump w/Clik-Tite Connectors for 2 pad use w/Mul-T-Pads (See GTA adaptor to use with K-Pads) | \$536.00 |

| Model | Circulating Warm Water Pads | w/Clik-Tite | w/Colder Style |
|-------|---|-------------|--------------------|
| LKP | Reusable, most durable vinyl 17" x 24" K-Pad Repairable – for single pad use | N/A | \$184.00 |
| DKP | Disposable (but washable) 18" X 24" K-Pad | N/A | 10/case @ \$204.00 |

Mul-T-Pads w/Clik-Tite Connectors are designed to permit two pads to run with one T/Pump
(Colder Style Connectors can only be used with one pad)

| | | | |
|--------|---|--------------------|--------------------|
| HD-MTP | Gaymar Mul-T-Pad 15" x 22" (1 or 2 pad use) Reusable, heavy duty, pliable polymer | \$62.00 | N/A |
| L-MTP | Gaymar Mul-T-Pad 18"x 26" (1 or 2 pad use) Lightweight polymer/non-woven fabric | 10/case @ \$134.00 | 10/case @ \$134.00 |
| I-MTP | Gaymar Mul-T-Pad 15" x 22" (1 or 2 pad use) Lightweight polymer/non-woven fabric | 20/case @ \$233.00 | 10/case @ \$133.00 |
| M-MTP | Gaymar Mul-T-Pad 13"x 18" (1 or 2 pad use) Lightweight polymer/non-woven fabric | 20/case @ \$230.00 | 10/case @ \$130.00 |
| S-MTP | Gaymar Mul-T-Pad 3" x 23" (1 or 2 pad use) Lightweight polymer/non-woven fabric | 10/case @ \$123.00 | N/A |

| Model | Accessories & Adaptors | Price/Case |
|-------|--|------------|
| TPS-1 | Four (4) caster mobile stand for T/Pump with storage basket | \$265.00 |
| BXK-2 | Removable Tubing Adaptor for KM 107 Pump or GTP-K Pump to allow use of Gaymar Mul-T-Pads | \$15.00 |
| GTAw | Gaymar Tubing Adaptor converts GTP Pump for use with K-Pads | \$28.00 |
| CT-F | Clik-Tite connector – female | \$3.00 |
| CS-F | Colder Style connector – female | \$3.00 |
| CT-M | Clik-Tite connector – male | \$3.00 |
| CS-M | Colder Style connector – male | \$3.00 |

Please Note: Prices do not include shipping.

