

Water Heating & Cooling Pumps and Pads

ThermoCare® offers a variety temperature controlled water pumps and pads for safe, consistent, heat therapy on the surgery table. Animals can be placed on the pads, which are attached to the heating/cooling circulator pump, during surgery to keep body temperature at the desired level and reduce risk of shock.

ThermoCare® offers the Gaymar® T-Pumps and Mul-T-Pads, Aquamatic® K-Modules and K-Pads, as well as various accessories and adaptors that can accommodate a range pumps and pads. Pads are either reusable and repairable, or disposable, and come in a variety of sizes. Optional slip-on fleece bags are machine washable for easy removal of blood, feces, or other fluids. ThermoCare® can also facilitate pump repairs; please call or email info@thermocare.com for more information.

T/Pump and Mul-T-Pads

Sandwich the patient in continuous warmth by connecting two pads to one pump

T/Pump Advanced Features

LED Display Panel and Selectable Temperature Settings

- Clear and intuitive, LED lights indicate current settings
- Classic T/Pump: Select from Cooling, Low Heat or High Heat
- Professional T/Pump: Select from 50°F, 95°F, 100°F or 107°F
- Solid light indicates temperature set point has been reached

Comfortable Handle

- Large grip area for easy transport of unit

Audible and Visual Alarm

- Alerts caregiver to flow, low water or internal issues

On/Standby Indicator

- Yellow light indicates unit has power
- Green light indicates unit is on

Translucent Reservoir

- Determine water level at a glance
- Reservoir has large capacity for less frequent filling
- Fill lines for both warming and cooling (water and ice)

Non-marking Rubber Feet

- Provide stability and shock absorption
- Feet easily realigned to fit T/Pump stand



Extra-Wide Leakproof Cap

- Simply fill with tap water for warming or tap water and ice for cooling and unit is ready for operation. Fill lines are indicated for both water and ice
- Tethered strap prevents cap loss




Secure Hose Attachments

- Permanently attached 10' hoses assure proper, secure connections
- Hoses can not be accidentally lost or discarded



Leakproof Connections

- Patented Clik-Tite® connectors snap into place quickly and easily
- Allows for two pads to be used simultaneously with one unit to treat multiple sites
- Non-restricted water flow



Accessory T/Pump Stand

- Provides ease of transporting to patient's room
- Storage in wire basket for pads and accessories



T/Flash Technology

- Slight temperature change stimulates patient's heat receptors
- Provides the "always there" feeling so the patient knows it's working

Safe, 3 Layer Safety System

- Dual temperature sensors
- Dual microprocessors
- Back-up mechanical thermostat to ensure safe operation

Safe, Effective Therapy

- When used in combination with our versatile Mul-T-Pads®, the T/Pump provides safe, effective localized therapy while easing patient pain and discomfort



Model	Water Heating Pumps	Price/Case†
KM-107	Aquamatic® K-Module 107 Pump for K-Pad w/10' Dual Tubing & Key (See BXK-2 adaptor to connect Gaymar® Mul-T-Pads to K-Module Pump)	\$614.00
GTP-K	Gaymar® T-Pump 650C w/Colder-style Connectors for K-Pads (See BXK-2 adaptor to use Gaymar® Mul-T-Pads)	\$486.00
GTP	Gaymar® T-Pump w/Clik-Tite Connectors for 2 pad use w/Mul-T-Pads (See GTA adaptor to use with K-Pads)	\$486.00

Model	Water Heating & Cooling Pumps (Heats & Cools, Tap Water)	Price/Case†
TP700C	Gaymar® T-Pump w/Colder-Style Connectors for K-Pads (See BXK-2 adaptor to use Gaymar® Mul-T-Pads)	\$536.00
TP700	Gaymar® T-Pump w/Clik-Tite Connectors for 2 pad use w/Mul-T-Pads (See GTA adaptor to use with K-Pads)	\$536.00

Model	Circulating Warm Water Pads	w/Clik-Tite†	w/Colder Style†
LKP	Reusable, most durable vinyl 17" x 24" K-Pad Repairable – for single pad use	N/A	\$184.00
DKP	Disposable (but washable) 18" X 24" K-Pad	N/A	10/case @ \$204.00

Mul-T-Pads w/Clik-Tite Connectors are designed to permit two pads to run with one T/Pump
(Colder Style Connectors can only be used with one pad)

HD-MTP	Gaymar® Mul-T-Pad 15" x 22" (1 or 2 pad use) Reusable, heavy duty, pliable polymer	\$62.00	N/A
L-MTP	Gaymar® Mul-T-Pad 18"x 26" (1 or 2 pad use) Lightweight polymer/non-woven fabric	10/case @ \$134.00	10/case @ \$134.00
I-MTP	Gaymar® Mul-T-Pad 15" x 22" (1 or 2 pad use) Lightweight polymer/non-woven fabric	20/case @ \$233.00	10/case @ \$133.00
M-MTP	Gaymar® Mul-T-Pad 13"x 18" (1 or 2 pad use) Lightweight polymer/non-woven fabric	20/case @ \$230.00	10/case @ \$130.00
S-MTP	Gaymar® Mul-T-Pad 3" x 23" (1 or 2 pad use) Lightweight polymer/non-woven fabric	10/case @ \$123.00	N/A

Model	Accessories & Adaptors	Price/Case†
TPS-1	Four (4) caster mobile stand for T/Pump with storage basket	\$265.00
BXK-2	Removable Tubing Adaptor for KM 107 Pump or GTP-K Pump to allow use of Gaymar® Mul-T-Pads	\$15.00
GTAw	Gaymar® Tubing Adaptor converts GTP Pump for use with K-Pads	\$28.00
CT-F	Clik-Tite connector – female	\$3.00
CS-F	Colder Style connector – female	\$3.00
CT-M	Clik-Tite connector – male	\$3.00
CS-M	Colder Style connector – male	\$3.00

†Please Note: prices do not include shipping

